

MOULDABLE RESTRAINTS

Adolescents - Infants - Pediatrics - New Born



RESTRAINING MATERIALS

SPLINTS are one of the materials which have to be made with most frequency in a hospital either by the nurses themselves or by the maintenance service. Either way, precious time is taken up during the working day in their construction.

This 2mil product is made in different sizes depending on the age of the patient, but all the products have the same basic concept:

- Flexibility to be moulded into either
 - INTRA-ARTERIAL
 - INTRAVENOUSpositions coupled with sufficient rigidity for maintaining the position.
- Includes a fastening device to keep the patient's hand away from his body
- Cushioning to avoid irritating the skin.
- No adhesives; fastening is by means of "VELCRO".
- Made entirely by high intensity ultrasonic welding to guarantee reusability.



Ref: FI-300
FI-305
FI-310
FI-311



FI-300

FI-305

FI-310

FI-311



ii The product is supplied with fastening straps; if required without fastening strap, put "sc" after the reference number !!



2mil artesanía sanitaria, s.l

✉: c/ de los Talleres nº 31 - 28430 ALPEDRETE (Madrid) SPAIN

☎: +34-916309208 FAX: +34-912793659 ER:+34-606957883

E-Mail: export@2mil.com

www.2mil.com

Distributor