

ACCUBOARD®

Disposable Neonatal, Infant
and Pediatric Limb-Boards



Neonatal & Infant
PRODUCTS

Comfort and Support for Baby

Ameritus AccuBoards feature a formable core that gently conforms to baby, allowing small limbs to be comfortably positioned and stabilized during I.V. procedures.

ACCUBOARD®

Disposable Neonatal, Infant and Pediatric Limb-Boards

Three Sizes to Fit Your Patients Two Styles to Fit your Preference

Features

- Ameritus **AccuBoards** feature an aluminum core that may be easily contoured to virtually any shape and position for maximum patient comfort and support
- Soft, lightweight fabric and durable foam pad provide a gentle touch to fragile skin
- Ameritus **AccuBoards** are also available with soft, foam straps for securing the board to the patient
- Ameritus **Limb-Boards** feature a soft, spun fabric covering, thick foam padding, and an inner core designed for additional comfort and support
- Latex Free and DEHP Free
- Made in the USA



Ameritus AccuBoard - Formable AB Series
Available with or without soft straps




Ameritus Limb-Board – FM Series

ORDERING INFORMATION

Catalog Number AccuBoard	Size Specifications	Catalog Number Limb-Board	Size Specifications
403AB	Neonatal 3/4" x 2-5/8" 1.91 cm x 6.67 cm	403-FM	Neonatal 1" X 3" 2.54 cm x 7.62 cm
403ABV	Neonatal (With Straps) 3/4" x 2-5/8" 1.91 cm x 6.67 cm	404-FM	Infant 1.5" x 4" 3.81 cm x 10.16 cm
404AB	Infant 1-5/8" x 3-7/8" 4.13 cm x 9.85 cm	405-FM	Pediatric 2" x 5" 5.08 cm x 12.70 cm
404ABV	Infant (With Straps) 1-5/8" x 3-7/8" 4.13 cm x 9.85 cm	Quantity (All Sizes): 48/Box	
405AB	Pediatric 1-7/8" x 4-7/8" 4.75 cm x 12.38 cm		
405ABV	Pediatric (With Straps) 1-7/8" x 4-7/8" 4.75 cm x 12.38 cm		



Manufactured for:
Kentec Medical, Inc.
17871 Fitch, Irvine, CA 92614
TEL 949.863.0810 FAX 949.833.9730
TOLL FREE (USA only) 800.825.5996
sales@ameritusmedical.com
www.kentecmedical.com



Neonatal & Infant
PRODUCTS

DISTRIBUTED BY:

