

# Face Masks

## Anatomical Face Masks



Mercury Medical's Face Masks are anatomically contoured to conform to natural facial profiles allowing a secure, low pressure seal without sharp edges that could damage the eyes.

### FEATURES:

- Completely disposable
- Clear material construction
- Anatomically designed
- Soft, pliable, air-filled cushion
- Removable color-coded headstrap fastener

### BENEFITS:

- Eliminates cross-contamination risks and associated reprocessing costs.
- Better patient visualization.
- Accommodates a wide range of facial features and reduces patient discomfort.
- Reduces wrinkles or creases and allows for a secure low pressure seal.
- Allows for use of traditional headstraps.

### ITEM NUMBER:

### DESCRIPTION:

10-50500	Adult Large Clear Anatomical Face Mask 30/case (Large)
10-50501	Adult Small Clear Anatomical Face Mask 30/case (Medium)
10-50502	Child Clear Anatomical Face Mask 20/case (Small)
10-50503	Pediatric Clear Anatomical Face Mask 20/case
10-50504*	Infant Clear Anatomical Face Mask 20/case
10-50505*	Neonatal Clear Anatomical Face Mask 20/case

\*Headstrap Fastener not available on these sizes.

## Round Silicone Masks

Features anatomically correct "average" contour for infant faces. Provides overall tight seal with internal cavity, minimizing dead space. Designed specifically for children under the age of ten (10).

10-50601	Size 0 10/pkg.
10-50602	Size 1 10/pkg.
10-50603	Size 2 10/pkg.

