



SunTech[®] CT50[™]

BP + VITALS

Vital Signs Monitor

Created with advanced monitoring features, the SunTech CT50 includes market-leading technology for all parameters as well as industry-standard connectivity to H.I.S. With the accuracy of BP averaging and the clarity of an intuitively designed, full-color, touch screen interface, this monitor has been created with both patients and clinicians in mind.

Standard Features

- 3 measurement modes (Monitor, Triage & STAT)
- Programmable monitoring intervals for BP & patient alarms
- 8" TFT Color Touch screen
- Masimo SET (preferred) or Nellcor SpO₂
- Covidien Filac 3000 Thermometry
- Ethernet LAN, USB x 4

Optional Accessories

- Thermal printer (internal)
- Bar code reader
- Mobile stand
- WiFi

