



DIMAR

Medical Devices

DimAir Line

The Helmet REVOLUTION

Intensive Care
Haematology
Infectivology
Home Care



ZIP HELMETS

Single patient Helmets for Non Invasive Mechanical Ventilation (NPPV), that can be used both in PSV-ASB, BIPAP and CPAP. Supplied with a pneumatic seal zip, they allow a quick and easy access to the patient. Comfortable and practical.

TECHNICAL FEATURES

- Opening Zip with pneumatic seal.
- Made of biocompatible plastic material.
- Additional fixing eyelets on the top of the helmet, for an alternative/additional connection
- Symmetrical right and left connectors for an easy and practical connection with the equipment placed on any side.
- Sealed access ports from 3 to 20 mm.
- Straps made of flexible soft material to reduce ulcerations
- Protective pads in medical approved antibacterial CoolMax® material
- 22M and 22F ISO connectors.
- Latex Free.



ADVANTAGES

- Quick access to the patient's face.
- In PSV-ASB modality, the fixing "eyelets" improve the inspiratory trigger reaction and the synchronism between the patient and the ventilator.
- Possibility to carry out cleaning, feeding and weaning operations without taking the helmet off.
- The gas inflow and outflow can be inverted or connected from the same side; in this way the patient can also be placed on his side.



NIMV 6R ZIP

Helmet for Non Invasive Mechanical Ventilation (NPPV), that can be used in PSV-ASB, BIPAP, and CPAP. 6 ISO connectors placed on the rigid ring for a greater Efficiency, Practicality and Comfort.

The connection to the rigid ring reduces the delay of the trigger reaction time <0,1".



COMFORT NEO-PED ZIP

CPAP Helmets for neonatal patients from 3 to 5 Kg. and paediatric patients from 5 to 15 Kg. Devoid of the rigid ring, they prevent the formation of ulceration to the Patient's neck. Fixing pad fabricated from

medical approved anti-bacterial Coolmax® material.



CPAP COMFORT ZIP

Helmet for CPAP. Devoid of the rigid ring, it prevents the formation of neck and shoulders ulcerations. The high comfort allows it to be used for prolonged periods both in hospital departments and in home care treatments.

NEONATAL	500/5447	3 -5 Kg.
PAEDIATRIC	500/5448	5 -15 Kg.

NIMV 6R ZIP	COMFORT CPAP ZIP	SIZE
600/5586	500/5579	XXS
600/5234	500/5441	XS
600/5235	500/5442	S
600/5236	500/5443	M
600/5237	500/5444	L
600/5238	500/5445	XL
600/5239	500/5446	XXL
600/5582	500/5568	3XL

- Single pack
- Multiple pack with 5 pieces



Possibility to realize different models, versions and packaging on customers' specific requirement.

DIMAR

Via G. Galilei, 6 41036 Medolla -MO- I
Tel: +39 0535 611336 Fax: +39 0535 611328
Mail: dimar.srl@tiscali.it

