

CURAPAX®

SINGLE-USE HEATING AND COOLING PADS

CURAPAX® CARE CONCEPT

This new and unique concept offers increased hygiene, increased accessibility and reduced treatment time. This brings hot and cold therapy directly into each and every hospital department, and makes treatment easier and more comfortable for both the patient and professional practitioners.



JUST ADD WATER



Compare	CURAPAX	Re-usable
Hygienic to use	√	-
Hygienic to store	√	-
Fast preparation	√	-
Use directly on skin	√	-
Space saving storage – hygienic dispenser box	√	-
Easily accessible	√	-
Integrated heatindicator	√	-

Specifications	
Rawmaterial	PE, PP, PA, SAP
Weight	22 grams
Size Curapax	250mm x 160mm
Size 50 unit package	110mm x 190mm x 285mm
Packaging	50 pcs./ box, 200 pcs. per carton
Latex and PVC free	-

Curapax single-use heat- and cold pad with integrated non-woven textile cover.

CURAPAX®

SINGLE-USE HEATING AND COOLING PADS

- Hygienic single-use concept which reduces cross contamination
- Timesaving preparation and disposal
- A safe solution for both patients and staff



FAST PREPARATION

The integrated one-way valve makes it easy to fill with water

MAXIMUM HYGIENE

Eliminate the risk of cross contamination with CURAPAX's disposable design

TIMESAVING

Increase efficiency with convenient access to pads, fast heating, and no cleanup. CURAPAX's non-woven exterior doesn't need to be wrapped in a towel

SAFE TO USE

Protect both patient and caregiver from burns and freezes with CURAPAX's nonwoven covering

LONGER TREATMENT TIME

Provide an extended treatment time due to CURAPAX's large capacity (550 ml gel) combined with its nonwoven insulating covering

FLEXIBILITY

It's easy to create customized applications such as neck/shoulder packs by connecting multiple packs via their unique corner 'eyes'.

INTEGRATED HEAT INDICATOR

Heating indicator is an additional indicative precaution against harmful use.

PRECAUTIONS:

If the intake valve is not completely empty a small amount of water/gel may leak. An overheated pack can burst – take care when removing the pack from the microwave/incubator. The temperature of the pack can continue to rise slightly after the heating has been terminated. As there is a heightened risk of using the product when the temperature reaches 45 degrees Celsius the heat indicator turns red. The heat indicator is solely an extra indicative precautionary measure against injurious use. The pad and the gel can be too warm for the desired usage - even if the indicator not has turned red. Please always handle the pad with caution at all times. A pad that is too warm or too cold can result in injuries on the skin and the underlying tissue. Take special care when using on patients with reduced sensitivity. People suffering from diabetes and people with poor circulation should use heat therapy only after consulting a physician. Be advised that the skin of children, the elderly and pregnant women can be especially sensitive why supervision is advised. The pack is ONLY for external use. If the pack is damaged it must be disposed of immediately.

RECOMMENDATION:

Each microwave performs differently. Therefore, we recommend you to heat CURAPAX at intervals of e.g. 20 seconds until you have gained a more precise sense of when the product reaches the desired temperature in your particular microwave.