



**New 26g**



Catetere intravenoso periferico  
con sistema di sicurezza passivo 26g  
**Neonatale**

*Intravenous Catheter with passive  
safety system 26g Neonatal*

DELTA SELFSAFE<sup>®</sup> **FEP**  
NEO DELTA SELFSAFE<sup>®</sup> **PUR**



New 26g



**Misura innovativa,**  
ideata con l'utilizzo dei migliori materiali  
per garantire il massimo comfort e precisione a operatori.  
Adatto a pazienti con vene di piccola dimensione.

*Innovative size, designed with the best materials  
to ensure the best comfort and precision for users.  
Suited for patients with smaller veins.*



- DELTA SELFSAFE - T
- NEO DELTA SELFSAFE - T



- DELTA SELFSAFE - 1
- NEO DELTA SELFSAFE - 1



- DELTA SELFSAFE - 2
- NEO DELTA SELFSAFE - 2

*Disponibile in 3 versioni*  
**Available in 3 Models**

FEP CODE	PUR CODE	DESCRIPTION	GAUGE	CATHETER LENGTH (mm) (inch)	Ø EXT. CATHETER (mm)	Ø INT. CATHETER (mm)	CAPACITY (ml/min)	SALES UNIT
3703222	3704022	Dual Entry I.V. Catheter with wings	26	19 +1	0.64	0.45	19	200
3706422	3708522	Single Entry I.V. Catheter without wings	26	19 +1	0.64	0.45	19	200
3706522	3708122	Single Entry I.V. Catheter with wings	26	19 +1	0.64	0.45	19	200

DELTA SELFSAFE® **FEP**

DELTA SELFSAFE® **PUR**