INFANT TRANSPORT MATTRESS



SAFE, PORTABLE HEATING



Our gel-based Infant Transport Mattress helps bring warm comfort to infants during transport in or between hospitals. Upon activation, the mattress reaches a maximum temperature and stays there, ensuring ultimate comfort and safety. Made with our patented Insta-Gel, our mattress is designed to comfortably reduce the risk of cold stress to infants whenever it's needed most, from inner-hospital transport to use in EMS vehicles.



BENEFITS

- **Safe and comfortable:** Gel-based for safe transport in or between hospitals
- Reduces cold stress: By keeping the infant warm during transport

FEATURES

- Single-use minimizes the risk of infection
- Portable heat source, no electricity required
- Delivers approximately 39.5°C/103°F of warm therapy when activated at room temperature and will not exceed 40.5°C/105°F
- · Easy activation disk

91650IM - Infant transport mattress is 9.65"x15.75



We have a variety of size and packaging options available to meet your needs.