


# with you all the way

safe | reliable | flexible



# Contents

About Sidam Group
Preparation
Connectors
Compounding
Administration







# From the past to the future through innovation

# The Sidam Group

In the city of the well-known humanist and philosopher Pico della Mirandola, in the Modena district, the biomedical district is currently the leader in the European market for the production of disposables and components for dialysis.

The Biomedical District of Mirandola includes over 100 companies specialized in disposable plastic and electro-medical products: for the most part they are Italian SMEs (85.6%) but there are also numerous multinationals (14.4%) both foreign ones (60%) and Italian (40%), which employ about 5,000 people and generate a turnover of almost 1 billion euros (including related industries).

Sidam was founded in January 1991 in the heart of the Mirandola district by Graziano Azzolini, an expert in the sector, who understands the need for international hospitals of excellence to have disposable medical devices, i.e., innovative and truly usable products.

Today Sidam is a company that covers all production phases, from R&D to after-sales activities. It currently works with 19 countries around the world and the thrust of exports continues to be the engine of growth.

In 2016 SIDAM acquired BTC Medical Europe, which was founded in 2006, with the focus on the production of devices for the oncology departments and medical devices for radio-pharmaceutical companies. BTC Medical Europe was fully incorporated by Sidam Group in March 2022.

4





Patent-protected, safety anti-disconnection system, for the connection to bags with cytostatic and chemotherapy drugs.

- Protected needle
- Safety anti-disconnection system
- Compatibility with all bag sizes
- Easy insertion through the secondary port of the bag

The Dilution operation, performed with non-protected needles, exposes the operator to the risk of accidental needlestick injuries.

The European Directive 2010/32/UE provides for the adoption of effective measures to protect workers against injuries caused by needles and other medical sharp devices in the hospital and healthcare sector

Code	ecoflac compatible	compatible with all bags	Antibac filter	silcone layer on tip	Vial Coupling system	Facilitated injection withdrawal maneuvers	Spike Protection Cap	Transparent
4033				•			•	•
4038	•	•	•				•	•
4047			•	•	20mm	•	•	•
4053			•	•	20mm	•	•	•
4090			•	•	13mm	•	•	•
4058			•	•	13mm	•	•	•
4062		•						•

caps



#### 4015 Smurfy

SMURFY guarantees a closed system at disconnection and it only activates when firmly connected to the female luer lock.

This prevents the escape of drug and vapours from the system. Once connected to the syringe, the Smurfy can rotate 360°.

Product specifications						
Flow rate by gravity at 91cm height	244 ml/min					
Functional activations	300					
Residual volume/dead space	0.09ml					







Supplied with protective cap

Supplied with protective cap

4010-S

The 4010-S is a closed circuit, high flow connector. With reinforced walls to avoid leakage with high pressures, the connector can rotate 360° and boasts an antidisconnection system and has an orange pigmentation to prevent UV damage.

COMPOUNDING COMPOUNDING

### Standard sets

Our standard systems eliminate the risk of drugs spreading during transport and protects the operator from contamination.



046-102 Standard Set



# Break off sets

A patent-protected system that guarantees the closure of the circuit until it is activated by the operator after connection to the administration set. This system eliminates the risks of drug spreading during transport and protects the operator.



Break off set



# Buddy sets


Patent-protected anti-disconnection system that allows irreversible connection to the administration set, preventing subsequent disconnection from it.



Buddy set

10 11

# Pumps/Gravity





These drugs infusion Sets are intended to be directly or indirectly connected to infusion pumps.

Code	Materials	Configuration	Technical Benefits	Qty/box
A-INCG01	PVC DEHP Free - Amber	1 way	back check valve One way	50
A-INCG02	PVC DEHP Free - Amber	2 ways	back check valve Two way	50
A-INCG02P	PVC DEHP Free - Amber	2 ways	Needle-free Two way	50
A-INCG04	PVC DEHP Free - Amber	4 ways	back check valve Four way	50
A-INCG04P	PVC DEHP Free - Amber	4 ways	Needle-free Four way	50
A-INCG06	PVC DEHP Free - Amber	6 ways	back check valve Six way	50
022-500	PVC DEHP Free - Amber	Onco Smurfy Set	Close Male Luer connector	50
046-600	PVC DEHP Free - Amber	IV Set	Y Injection Site	200
046-750	PE - Amber	IV Set	Compatibility with all drugs	200

All sets are also available with transparent tubes

12 13



safe | reliable | flexible

Via Statale Sud, 169 41037 - Mirandola (MO) ITALY

www.sidamgroup.com

Design: leecallaghan.com | Printed by: Tipolito Salvioli